

Three Line Acupuncture for Bell's Palsy

58 cases of Bell's Palsy patients were treated with this method during the last 7 years. Among all the cases, 52 cases were cured. The total effective rate is 100%.

1. Location of the three lines

	Starting points	Ending points
Line one	Taiyang(extra)-----	Tinggong (SI19)
Line two	Yingxiang(LI20)-----	Xiaguan(ST7)
Line three	Dicang(ST4)-----	Jiache(ST6)

2. Method

Needle the starting, middle and the ending points of each line. Nine points all together. Give even manipulation (not strong stimulation), and retain the needle for 20 minutes each treatment. Stimulate the needles every 5 minutes during the retention. The treatment is given every day and lasts for 2 courses. One course is 7 days.